Embedding of *Oxyrrhis marina* for TEM (by Erhard Rhiel)

Day 1

- harvest the cells from 800–1000 ml of culture medium by centrifugation (we use an Eppendorf 5810R refrigerating centrifuge equipped with an A-4-62 swinging bucket rotor, and centrifuge for 10 min at 100×g)
- fix the cells for 2 h at room temperature using either 4 % (v/v) glutaraldehyde in 0.1 M Na–K phosphate buffer (pH 7) containing 0.05 % (w/v) ruthenium red or 4 % (v/v) glutaraldehyde in f/2 medium containing 20 mM Na–K phosphate buffer (pH 7) and 0.05 % (w/v) ruthenium red
- wash the cells four/five times with 0.1 M Na–K phosphate buffer (pH 7); centrifuge them in between (we use a rather old table top centrifuge equipped with a swinging bucket rotor, and centrifuge for 5 min at 1000 rpm)
- perform a post-fixation step for 1.5–2 h in 1 % (v/v) OsO₄
- wash the cells four/five times with 0.1 M Na–K phosphate buffer (pH 7); centrifuge them in between
- dehydrate the cells in a graded ethanol series. Note: you can add 2 % (w/v) uranyl acetate during the 50 % ethanol step
- 30 min., 30 % EtOH
- 60 min., 50 % EtOH with 2% Uranylacetat
- 30 min., 70 % EtOH
- over night, 70 % EtOH

Day 2

- 30 min., 80 % EtOH
- 30 min., 90 % EtOH
- 30 min., 100 % EtOH
- 1 h, 100 % EtOH
- 1 h, EtOH:Spurr (2:1)
- 2 h, EtOH:Spurr (1:1)
- 2 h, EtOH:Spurr (1:2)
- 2 h, Spurr resin (pure)
- Incubate over night in fresh Spurr resin (pure) at room temperature and under the hood

Day 3 (Embedding of the samples in epoxy resin according to Spurr)

- centrifuge for 10 min and at 1400 rpm
- resuspend in fresh Spurr in small cups
- spin down the cells to the bottom by centrifugation
- incubate 2h at 40 °C followed by 70 °C over night

Literature:

Spurr AR, 1969: A low-viscosity epoxy resin embedding medium for electron microscopy. J Ultrastruct Res 26:31–43

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